

## Product datasheet

# Recombinant Mouse CD134 / OX40L receptor protein (His tag) ab215612

### Description

---

<b>Product name</b>	Recombinant Mouse CD134 / OX40L receptor protein (His tag)
<b>Purity</b>	> 95 % SDS-PAGE.
<b>Endotoxin level</b>	< 1.000 Eu/μg
<b>Expression system</b>	Mammalian
<b>Accession</b>	<b><u>P47741</u></b>
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Mouse
<b>Sequence</b>	VTARRLNVCVKHTYPSGHKCCRECQPGHGMVSRCDHTRD TLCHPCETGFYN EAVNYDTCKQCTQCNRHRSSELKQNCPTQDTVCRCRP GTQPRQDSGYKL GVDCVPCPPGHFSPGNNQACKPWTNCTLSGKQTRHPAS DSLDAVCEDRSL LATLLWETQRPTFRPTTVQSTTVWPRTSELPSPTLVTP GPHHHHHH
<b>Predicted molecular weight</b>	22 kDa including tags
<b>Amino acids</b>	20 to 211
<b>Tags</b>	His tag C-Terminus

### Specifications

---

Our **Abpromise guarantee** covers the use of **ab215612** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** SDS-PAGE

**Form** Lyophilized

### Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Store at -20°C.

pH: 7.40

Constituent: 100% PBS

**Reconstitution** Reconstitute with ddH<sub>2</sub>O to a concentration no less than 100 µg/mL. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

## General Info

---

**Function** Receptor for TNFSF4/OX40L/GP34.

**Sequence similarities** Contains 4 TNFR-Cys repeats.

**Cellular localization** Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors