

## Product datasheet

# Recombinant Mouse CD134 / OX40L receptor protein (Tagged) ab271643

1 Image

### Description

---

<b>Product name</b>	Recombinant Mouse CD134 / OX40L receptor protein (Tagged)
<b>Purity</b>	>= 90 % SDS-PAGE.
<b>Expression system</b>	HEK 293 cells
<b>Accession</b>	<b><u>P47741</u></b>
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Mouse
<b>Sequence</b>	VTARRLNCVKHTYPSGHKCCRECCQPGHGMVSRCDHTRD TLCHPCETGFYN EAVNYDTCKQCTQCNHRSGSELKQCNCTPTQDTVCRCRP GTQPRQDSGYKL GVDCVPCPPGHFSPGNNQACKPWTNCTLSGKQTRHPAS DSLDAVCEDRSL LATLLWETQRPTFRPTTVQSTTVWPRTSELPSPTLVTP GP
<b>Predicted molecular weight</b>	75 kDa
<b>Molecular weight information</b>	This protein runs at a higher MW by SDS-PAGE due to glycosylation.
<b>Amino acids</b>	20 to 211
<b>Tags</b>	Avi tag C-Terminus , Fc tag C-Terminus
<b>Additional sequence information</b>	Extracellular domain.

### Specifications

---

Our **Abpromise guarantee** covers the use of **ab271643** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** SDS-PAGE

**Form** Liquid

## Preparation and Storage

---

### Stability and Storage

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 0.13% Sodium phosphate, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

## General Info

---

### Function

Receptor for TNFSF4/OX40L/GP34.

### Sequence similarities

Contains 4 TNFR-Cys repeats.


### Cellular localization

Membrane.

---

## Images

---



SDS-PAGE analysis of 4 µg ab271643.

This protein runs at a higher MW by SDS-PAGE due to glycosylation.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors