

Product datasheet

Recombinant Mouse CD5L/CT-2 protein (His tag)
ab208467

1 Image

Description	
Product name	Recombinant Mouse CD5L/CT-2 protein (His tag)
Purity	> 95 % SDS-PAGE. Affinity purified
Endotoxin level	< 1.000 Eu/µg
Expression system	Insect cells
Accession	<u>Q9QWK4</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Mouse
Sequence	ESPTKVQLVG GAHRCEGRVE VEHNGQWGTV CDDGWDRRDV AVVCRELNCG AVIQTPRGAS YQPPASEQRV LIQGVDCNGT EDTLAQCELN YDVFDCSHEE DAGAQCENPD SDLLFIPEDV RLVDGPGHCQ GRVEVLHQSQ WSTVCKAGWN LQVSKVVC RQ LGCGRALLTY GSCNKSTQGK GPIWMGKMSC SGQEANLRSC LLSRLENNCT HGEDTWMECE DPFELKLVGG DTPCSGRLEV LHKGSWGSVC DDNWGEKEDQ VVCKQLGCGK SLHPSPKTRK IYGPAGARIW LDDVNCSGKE QSLEFCRHRL WGYHDC THKE DVEVICTDFD VLEHHHHHH
Predicted molecular weight	38 kDa including tags
Amino acids	22 to 352
Tags	His tag C-Terminus
Additional sequence information	NCBI Accession No.: AAH94459. SDS–PAGE under reducing conditions the protein migrates to 40-57 KDa.

Specifications

Our **Abpromise guarantee** covers the use of **ab208467** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid
Additional notes	Previously labelled as CD5L.


Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.40 Constituents: 90% PBS, 10% Glycerol (glycerin, glycerine)
------------------------------	---

General Info

Function	May play a role in the regulation of the immune system. Seems to play a role as an inhibitor of apoptosis.
Tissue specificity	Expressed in spleen, lymph node, thymus, bone marrow, and fetal liver, but not in non-lymphoid tissues.
Sequence similarities	Contains 3 SRCR domains.
Cellular localization	Secreted.

Images



15% SDS-PAGE analysis of ab208467 (3 µg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors