

Recombinant Mouse CD62L protein ab205901

1 Image

Description

Product name	Recombinant Mouse CD62L protein		
Purity	> 95 % SDS-PAGE.		
Expression system	Escherichia coli		
Accession	<u>P18337</u>		
Protein length	Protein fragment		
Animal free	No		
Nature	Recombinant		
Species	Mouse		
Sequence	QCQYVVQCEPLEAPELGTMDCIHPLGNFSFQSKCAFNCSEGRELLGTAETQCGASGNWSSPEPICQVTVQCEPLEAPELGTMDCIHPLGNFSFQSKCAFNCSEGRELLGTAETQCGASGNWSSPEPICQETNRSFSKIKEGDYN		
Predicted molecular weight	20 kDa including tags		
Amino acids	190 to 332		

Specifications

Our **Abpromise guarantee** covers the use of **ab205901** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
	Sandwich ELISA
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.
	pH: 8.00
	Constituents: 0.29% Sodium chloride, 0.01% EDTA, 0.61% Tris, 1% Glycine, 8.7% DL-Arginine,

General Info

Function	Cell surface adhesion protein. Mediates the adherence of lymphocytes to endothelial cells of high endothelial venules in peripheral lymph nodes. Promotes initial tethering and rolling of leukocytes in endothelia.
Tissue specificity	Expressed in B cell lines and T lymphocytes.
Sequence similarities	Belongs to the selectin/LECAM family. Contains 1 C-type lectin domain. Contains 1 EGF-like domain. Contains 2 Sushi (CCP/SCR) domains.
Cellular localization	Membrane.

Images



SDS-PAGE analysis of ab205901 (1 µg) stained with colloidal coomassie blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors