

Product datasheet

Recombinant Mouse GST3 / GST pi protein ab219251

1 Image

Description

Product name	Recombinant Mouse GST3 / GST pi protein
Purity	> 95 % SDS-PAGE. ab219251 was purified by using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>P46425</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Mouse
Sequence	MPPYTIVYFPSPGRCEAMRMLLADQQQSWKEEVVTIDTW MQGLLKPTCLY GQLPKFEDGDLTLYQSNAILRHLGRSLGLYGKNQREAAQV DMVNDGVEDL RGKYGTMIYRNYENGKNDYVKALPGHLKPFETLLSQNQGG KAFVGDQIS FADYNLLDLLLIHQVLAPGCLDNFPLLSAYVARLSARPKIK AFLSSPEHV NRPINGNGKQ
Predicted molecular weight	24 kDa
Amino acids	1 to 210

Specifications

Our **Abpromise guarantee** covers the use of **ab219251** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------------------------	--


pH: 7.40

Constituents: 20% Glycerol (glycerin, glycerine), PBS

General Info

Function	Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.
Sequence similarities	Belongs to the GST superfamily. Pi family. Contains 1 GST C-terminal domain. Contains 1 GST N-terminal domain.

Images



15% SDS-PAGE analysis of ab219251 (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors