

Recombinant Mouse Thyroxine Binding Globulin protein (His tag) ab276810

1 Image

Description

Product name	Recombinant Mouse Thyroxine Binding Globulin protein (His tag)	
Purity	> 97 % SDS-PAGE.	
Endotoxin level	< 1.000 Eu/µg	
Expression system	HEK 293 cells	
Accession	<u>P61939</u>	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Mouse	
Sequence	MSVFFYLFVL VFGLQATIHC APHNSSEGKV TTCHLPQQNA TLYKMPSINA DFAFSLYRRL SVENPDLNI FFSPVSISVA LAMLSFGSGS STQTQILEVL GFNLTDTPVT ELQQGFQHLL CSLNFPKNE LELQMGNAVF IGQQLKPLAK FLDDVKTLYE TEVFSTDFSN VSAAQHKINS YVEKQTKGK VGLIQGLKL NIIMILVNYI HFRAQWANPF RVSKTEESSN FSVDKSTTVQ VPMMHQLEQ YYHYVDMELN CTVLQMDYSE NALALFVLPK EGHMEWVEAA MSSKTLKKWN YLLQKGWVE LFVPKFSSA TYDLGSTLQK MGMRDAFAES ADFPGITEDS GLKLSYAFHK AVLHIGEEG TKEGASPEVG SLDQQEVPPL HPVIRLDRAF LLMILEKRTR SVLFLGKLVN PTKQ	
Predicted molecular weight	46 kDa including tags	
Molecular weight information	The secreted recombinant mouse SERPINA7 comprises 409 amino acids and has a predicted molecular mass of 46.2 kDa. As a result of glycosylation, it migrates as an approximately 55-60 kDa band in SDS-PAGE under reducing conditions.	
Amino acids	1 to 418	
Tags	His tag C-Terminus	
Additional sequence information	Predicted N-terminus: Ala 21.	

Specifications

Our **Abpromise guarantee** covers the use of **ab276810** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Lyophilized

Preparation and Storage

Stability and Storage Shipped at Room Temperature. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.


pH: 7.40

Constituent: 100% PBS

Reconstitution This information is lot specific. Please contact our technical Support team for details.

General Info

Images



SDS-PAGE analysis of ab276810

SDS-PAGE - Recombinant Mouse Thyroxine
Binding Globulin protein (His tag) (ab276810)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors