

Recombinant Mouse TIE2 protein ab206798

1 Image

Description	
Product name	Recombinant Mouse TIE2 protein
Purity	> 95 % SDS-PAGE. ab206798 was expressed in insect cell and affinity purified.
Endotoxin level	< 1.000 Eu/µg
Expression system	Insect cells
Accession	<u>Q02858</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Mouse
Sequence	AMDILINSLPLVSDAETSLTCIASGWHPEPITIGRDFEAL MNQHQPDL EVTQDVTREWAKKVWVKREKASKINGAYFCEGRVRGQAI RIRTMKMRQQA SFLPATLTMTVDRGDNVNISFKKVLKEEDAVMYKNGSFIHS VPRHEVPD ILEVHLPHAQPQDAGVYSARYIGGNLFTSAFTRLIVRRCEA QKWGPDCSR PCTTCKNNGVCHEDTGECICPPGFMGRTCEKACEPHTFG RTCKERCSGPE GCKSYVFCLPDPYGCSCATGWRGLQCNEACPSGYYGPD CKLRCHCTNEEI CDRFQGCLCSQGWWQGLQCEKEGRPRMTPQIEDLPDIE VNSGKFNPICKA SGWPLPTSEEMTLVKPDGTVLQPNDFNITDRFSVAIFTVN RVLPPDSGVW VCSVNTVAGMVEKPFNISVKVLPEPLHAPNVIDTGHNFAII NISSEPYFG DGPIKSKKLFYKPVNQAWKYIEVTNEIFTLNILEPRTDYELC VQLARPGE GGEGHPGPVRRFTTASIGLPPPRGLSLLPKSQTALNLTWQ PIFTNSEDEF YVEVERRSLQTTSDQQNIKVPGNLTSVLLSNLVPREQYTV RARVNTKAQG

EWSEELRAWTLSDILPPQPENIKISNITDSTAMVSWTIVDGY
 SSSIIIR
 YKVQGKNEDQHIDVKIKNATVTQYQLKGLEPETTYHVDIFA
 ENNIGSSNP AFSHELRTLPHSPASADLGGGKMLHHHHHH

Predicted molecular weight	82 kDa including tags
Amino acids	23 to 746
Tags	His tag C-Terminus
Additional sequence information	Extracellular domain. NP_001277478.

Specifications

Our **Abpromise guarantee** covers the use of **ab206798** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.40 Constituents: 90% PBS, 10% Glycerol (glycerin, glycerine)
------------------------------	--

General Info

Function	Tyrosine-protein kinase that acts as cell-surface receptor for ANGPT1, ANGPT2 and ANGPT4 and regulates angiogenesis, endothelial cell survival, proliferation, migration, adhesion and cell spreading, reorganization of the actin cytoskeleton, but also maintenance of vascular quiescence. Has anti-inflammatory effects by preventing the leakage of proinflammatory plasma proteins and leukocytes from blood vessels. Required for normal angiogenesis and heart development during embryogenesis. Required for post-natal hematopoiesis. After birth, activates or inhibits angiogenesis, depending on the context. Inhibits angiogenesis and promotes vascular stability in quiescent vessels, where endothelial cells have tight contacts. In quiescent vessels, ANGPT1 oligomers recruit TEK to cell-cell contacts, forming complexes with TEK molecules from adjoining cells, and this leads to preferential activation of phosphatidylinositol 3-kinase and the AKT1 signaling cascades. In migrating endothelial cells that lack cell-cell adhesions, ANGPT1 recruits TEK to contacts with the extracellular matrix, leading to the formation of focal adhesion complexes, activation of PTK2/FAK and of the downstream kinases MAPK1/ERK2 and MAPK3/ERK1, and ultimately to the stimulation of sprouting angiogenesis. ANGPT1 signaling triggers receptor dimerization and autophosphorylation at specific tyrosine residues that then serve as binding sites for scaffold proteins and effectors. Signaling is modulated by ANGPT2 that has lower affinity for TEK, can promote TEK autophosphorylation in the absence of ANGPT1, but inhibits ANGPT1-mediated signaling by competing for the same binding site. Signaling is also modulated by formation of heterodimers with TIE1, and by proteolytic processing that gives rise to a soluble TEK extracellular domain. The soluble extracellular domain modulates signaling by functioning as decoy receptor for angiopoietins. TEK phosphorylates DOK2, GRB7, GRB14, PIK3R1; SHC1 and TIE1.
-----------------	--

Tissue specificity	Detected in umbilical vein endothelial cells. Proteolytic processing gives rise to a soluble extracellular domain that is detected in blood plasma (at protein level). Predominantly expressed in endothelial cells and their progenitors, the angioblasts. Has been directly found in placenta and lung, with a lower level in umbilical vein endothelial cells, brain and kidney.
Involvement in disease	Dominantly inherited venous malformations May play a role in a range of diseases with a vascular component, including neovascularization of tumors, psoriasis and inflammation.
Sequence similarities	Belongs to the protein kinase superfamily. Tyr protein kinase family. Tie subfamily. Contains 3 EGF-like domains. Contains 3 fibronectin type-III domains. Contains 2 Ig-like C2-type (immunoglobulin-like) domains. Contains 1 protein kinase domain.
Domain	The soluble extracellular domain is functionally active in angiopoietin binding and can modulate the activity of the membrane-bound form by competing for angiopoietins.
Post-translational modifications	Proteolytic processing leads to the shedding of the extracellular domain (soluble TIE-2 alias sTIE-2). Autophosphorylated on tyrosine residues in response to ligand binding. Autophosphorylation occurs in trans, i.e. one subunit of the dimeric receptor phosphorylates tyrosine residues on the other subunit. Autophosphorylation occurs in a sequential manner, where Tyr-992 in the kinase activation loop is phosphorylated first, followed by autophosphorylation at Tyr-1108 and at additional tyrosine residues. ANGPT1-induced phosphorylation is impaired during hypoxia, due to increased expression of ANGPT2. Phosphorylation is important for interaction with GRB14, PIK3R1 and PTPN11. Phosphorylation at Tyr-1102 is important for interaction with SHC1, GRB2 and GRB7. Phosphorylation at Tyr-1108 is important for interaction with DOK2 and for coupling to downstream signal transduction pathways in endothelial cells. Dephosphorylated by PTPRB. Ubiquitinated. The phosphorylated receptor is ubiquitinated and internalized, leading to its degradation.
Cellular localization	Cell membrane. Cell junction. Cell junction, focal adhesion. Cytoplasm, cytoskeleton. Secreted. Recruited to cell-cell contacts in quiescent endothelial cells. Colocalizes with the actin cytoskeleton and at actin stress fibers during cell spreading. Recruited to the lower surface of migrating cells, especially the rear end of the cell. Proteolytic processing gives rise to a soluble extracellular domain that is secreted.

Images



15% SDS-PAGE analysis of ab206798 (3 µg).

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors