

Recombinant NCR1 protein ab51221

[1 References](#) [1 Image](#)

Description

Product name	Recombinant NCR1 protein
Purity	> 95 % SDS-PAGE. ab51221 was purified by FPLC gel-filtration chromatography, after refolding of the isolated inclusion bodies in a redox buffer.
Expression system	Escherichia coli
Protein length	Protein fragment
Animal free	No
Nature	Recombinant

Specifications

Our **Abpromise guarantee** covers the use of **ab51221** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage


Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 7.40 Constituents: PBS, 0.0292% EDTA
------------------------------	--

General Info

Function	Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.
Tissue specificity	Selectively expressed by both resting and activated NK cells.
Sequence similarities	Belongs to the natural cytotoxicity receptor (NCR) family. Contains 2 Ig-like (immunoglobulin-like) domains.
Post-translational	N-glycosylated.

modifications	O-glycosylated.
Cellular localization	Cell membrane.

Images



SDS analysis of recombinant extracellular Ig like domain of human NCR (ab51221). Gel: 14%.

SDS-PAGE - Recombinant NCR1 protein (ab51221)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors