

Recombinant Nucleoprotein (His tag) ab224846

1 Image

Description

Product name Recombinant Nucleoprotein (His tag)
Purity > 90 % SDS-PAGE.

Expression system Escherichia coli
Accession **P13699**
Protein length Full length protein
Animal free No
Nature Recombinant
Species Lassa virus

Sequence MSASKEIKSFLWTQSLRRELSGYCSNIKLVVVKDAQALLH
GLDFSEVSNV
QRLMRKERRDDNDLKRLRDLNQAVNNLVELKSTQQKSIL
RVGTLTSDDLL
ILAADLEKLKSKVIRTERPLSAGVYMGNLSSQQLDQRRALL
NMIGMSGGN
QGARAGRDGVVRVWDVKNAELLNNQFGTMPSLTLACLT
KQGQVDLNDVQ
ALTDLGLIYAKYPNTSDLDRLTQSHPIILNMIDTKKSSLNISG
YNFSLGA
AVKAGACMLDGGNMLETIKVSPQTMGILKSILKVKKALG
MFISDTPGER
NPYENILYKICLSGDGWPIASRTSITGRAWENTVVDLESD
GKPQKADSN
NSSKSLQSAGFTAGLTYSQMLTKDAMLQLDPNAKTWMD
IEGRPEDPVEI
ALYQPSSGCYIHFFREPTDLKQFKQDAKYSHGIDVTLFAT
QPGLTSAVI
DALPRNMVITCQGSDDIRKLLESQGRKDIKLIDIALSKTDSR
KYENAVWD
QYKDLCHMHTGVVVEKKKRGKKEITPHCALMDCIMFDA
AVSGGLNTSVL RAVLPRDMVFRTSTPRVVL

Predicted molecular weight 79 kDa including tags
Amino acids 1 to 569

Tags	His tag N-Terminus
Additional sequence information	Lassa virus (strain Mouse/Sierra Leone/Josiah/1976) (LASV) protein with 6xHis-SUMO tag at the N-terminus.

Specifications

Our **Abpromise guarantee** covers the use of **ab224846** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid


Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.2

Constituents: 50% Glycerol (glycerin, glycerine), Tris buffer

Images



Discontinuous SDS-PAGE (Tris-Glycine gel)(reduced) analysis of ab224846 with 5% enrichment gel and 15% separation gel.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors