

Product datasheet

Recombinant Pneumolysin protein (Tagged) ab236193

1 Image

Description

Product name	Recombinant Pneumolysin protein (Tagged)	
Purity	> 90 % SDS-PAGE.	
Expression system	Escherichia coli	
Accession	<u>P0C2J9</u>	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Streptococcus	
Sequence	ANKAVNDFILAMNYDKKKLLTHQGESIENRFIKEGNQLPDE FVVIERKKR SLSTNTSDISVTATNDSRLYPGALLVVDETLLENNPTLLAV DRAPMTYSI DLPGLASSDSFLQVEDPSNSSVRGAVNDLLAKWHQDYG QVNNVPARMQYE KITAHSMQLKVKFGSDFEKTGNSLDIDFNSVHSGEKQIQI VNFKQIYT VSVDVAVKNPGDVFQDTVTVEDLKQRGISAERPLVYISSVA YGRQVYLKLE TTSKSDEVEAAFEALIKGVKVAPQTEWKQILDNTEVKAVIL GGD PSSGAR VVTGKVDMMVEDLIQEGSRFTADHPGLPISYTTSLFRDNVV ATFQNSTDYV ETKVTAYRNGDLLLDHSGAYVAQYYITWNELSYDHQKEV LTPKAWDRNG QDLTAHFTTISIPLKGNVRNLSVKIRECTGLAWEWVRTVYE KTDLPLVRKR TISIWGTTLYPQVEDKVEND	
Predicted molecular weight	69 kDa including tags	
Amino acids	2 to 471	
Tags	His tag N-Terminus	
Additional sequence information	Streptococcus pneumoniae serotype 4. N-terminal 6xHis-SUMO-tagged.	

Specifications

Specifications

Our **Abpromise guarantee** covers the use of **ab236193** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle.


pH: 7.2

Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)

General Info

Cellular localization Plasma membrane

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) analysis with 5% enrichment gel and 15% separation gel of ab236193.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors