

Product datasheet

Recombinant *S. mansoni* Major egg antigen protein (His tag) ab224849

1 Image

Description

Product name	Recombinant <i>S. mansoni</i> Major egg antigen protein (His tag)	
Purity	> 90 % SDS-PAGE.	
Expression system	Escherichia coli	
Accession	<u>P12812</u>	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	<i>Schistosoma mansoni</i>	
Sequence	MSGGKQHNAVSIIPVNREQRSFEKQRRDLLTGLEHGGGAH RGNSIAPYTED WPSTVDNWDSSWKRWDDDMRRLRRGMFALLPLDTFSI GILENPFALMHQ MDRQIQDIRERMGSLDVPSTGSVNDFLKDAYEVEGEDGKV HFKVRFDAQGF APQDINVTSSENRVTVHAKKETTTDGRKCSREFCRMVQL PKSIDDSQLKC RMTDDGVLMLLEAPVKVDQNSLTLNESGQVAVRPKSDN QIKAVPASQALV AKGVHGLSYDDGSGGKRLHVEVAVDPVYKPEDLFVNV DSNRVVVSGRHH KQKSDQHGRSSSFAEFSQS YAIPETVDPLSVSAQVVGNT LVLEAPLEKQH AITH	
Predicted molecular weight	43 kDa including tags	
Amino acids	1 to 354	
Tags	His tag N-Terminus	
Description	Recombinant <i>S. mansoni</i> Major egg antigen protein (His tag)	

Specifications

Our **Abpromise guarantee** covers the use of **ab224849** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid


Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.2

Constituents: 50% Glycerol (glycerin, glycerine), Tris buffer

Images



Discontinuous SDS-PAGE (Tris-Glycine gel) (reduced) analysis of ab224849 with 5% enrichment gel and 15% separation gel.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors