



Product datasheet

20x TBS-T with Tween 20 ab64204

★★★★★ [1 Abreviews](#) [1 References](#) [1 Image](#)

Overview

Product name	20x TBS-T with Tween 20
Tested applications	Suitable for: IHC-P
General notes	<p>Tris Buffered Saline (TBS) wash buffer with Tween 20 for use during immunolabeling techniques such as IHC.</p> <p>The salts in 20x solution are at a very high concentration. It has been found that at temperatures below 44 degrees F (6 degrees C), the salts fall out of solution. If this occurs during your shipment conditions, please stir the product to get the salts back into solution.</p>

Properties


Form	Liquid
Storage instructions	Store at room temperature.
Storage buffer	<p>pH: 7.2</p> <p>Constituents: 99% Tris buffered saline (x20), 1% Tween</p> <p>20x Concentrated</p>
Relevance	<p>TBS is a solution of Tris Buffered Saline for use with immunolabeling techniques as a wash buffer. Tween 20 is sometimes added to promote more effective washings, resulting in decreased non-specific background staining. Tween 20 is also known to promote reagent spreading in automated staining methodologies.</p>

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab64204 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent dilution.

Images



ab64204 20X TBS with Tween 20 used along with an anti-CD20 antibody on Human Tonsil

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - 20x TBS-T with Tween 20 (ab64204)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors