

Product datasheet

BSA protein standards ab270701

Overview

Product name	BSA protein standards
General notes	<p>This product is native (Bovine) and is unpurified.</p> <p>Precisely formulated at 2mg/mL and are designed for serial dilution to construct standard curves.</p> <p>10 x 1 ml at 2 mg/mL.</p> <p>This product is manufactured by Expedeon, an Abcam company. It was previously called BSA Standards for quantification reagents 10 x 1ml (2mg/ml). Expedeon product code BS001 is the same as the 10 unit size of this product.</p>

Properties

Form	Liquid
Storage instructions	Shipped at Room Temperature. Store at -20°C long term.
Storage buffer	<p>pH: 7</p> <p>Constituent: 0.2% BSA</p>

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors