

## Product datasheet

### Cell Lysis Buffer ab288311

#### Overview

---

<b>Product name</b>	Cell Lysis Buffer
<b>General notes</b>	<p>This product is manufactured by BioVision, an Abcam company and was previously called 1067 Cell Lysis Buffer. 1067-100 is the same size as the 100 ml size of ab288311.</p> <p>The Cell Lysis Buffer (ab288311) is optimized for lysis of mammalian cells and tissues for caspase activity assays. The buffer is intended to use as a supplemental reagent for researchers who need to lyse large amount of cells or tissues.</p>

#### Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 7.2 Constituents: 0.02% (R*,R*)-1,4-Dimercaptobutan-2,3-diol, 1% Alcohol Ethoxylate

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors