

Cresyl Violet Stain Solution (1%) ab246817

[1 References](#) [1 Image](#)


Overview

Product name	Cresyl Violet Stain Solution (1%)
General notes	<p>Cresyl Violet Stain Solution (1%) ab246817 is designed for staining Nissl substance in neurons on formalin fixed, paraffin-embedded tissue. Nissl granules are purple, nuclei of neuroglia and endothelial cells are slightly bluer than Nissl granules.</p> <p>Nissl Substance: Violet.</p> <p>Neuroglia Nuclei: Violet to Dark Blue.</p> <p>Control Tissue: Cerebral Cortex</p> <p>Procedure:</p> <p>Working Solutions: If preferred, reagent may be diluted up to 1:10 with Deionized water just before use.</p> <ol style="list-style-type: none"> 1. Deparaffinize sections and hydrate to distilled water. 2. Apply Cresyl Violet Acetate solution (or working solution) to tissue for 3-5 minutes 3. Rinse quickly in 1 change of distilled water. 4. Dehydrate quickly in absolute alcohol. (Alcohol may remove stain from tissue over time). 5. Clear in Xylene or Xylene substitute and mount in synthetic resin. <p>Other products for staining tissue sections</p> <p>Find more kits and reagents in the special stains guide, or products for antigen retrieval, blocking, signal amplification, visualization, counterstaining, and mounting in the IHC kits and reagents guide.</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	<p>pH: 3.00</p> <p>Constituents: 0.2% Acetic acid, 1% 3,7-diamino-12H-benzo[b]phenoxazin-12-ium acetate</p>

Images



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Cresyl Violet Stain Solution (1%) (ab246817)

ab246817 Cresyl Violet Stain Solution (1%) staining formalin-fixed-paraffin embedded human brain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors