

## Product datasheet

### Normal Hamster Serum (sterile) ab7483

#### Overview

<b>Product name</b>	Normal Hamster Serum (sterile)
<b>Host species</b>	Syrian hamster
<b>Tested applications</b>	<b>Suitable for:</b> Blocking, IHC-P, IHC-Fr
<b>General notes</b>	Total Protein Concentration : normal Strain: Golden Syrian Sex: Mixed sex Raised in: Syrian Hamster Purity: Whole antiserum

#### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

#### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab7483 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Blocking		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors