

Product datasheet

Perchloric Acid (70%) ab291263

Overview

Product name	Perchloric Acid (70%)
General notes	This product is designed for use with Hydroxyproline Assay Kit ab222941 / K226, in combination with ab291264 (M1526) to produce the optional Perchloric acid/Isopropanol Solution that may be used instead of the Developer component.

Properties

Form	Liquid
Storage instructions	Shipped at Room Temperature. Store at Room Temperature. Store In the Dark.
Storage buffer	Constituent: 70% Perchloric acid

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors