

Protein A/G Magnetic Beads ab286842

Overview

Product name	Protein A/G Magnetic Beads
Tested applications	Suitable for: IP
General notes	<p>This product is manufactured by BioVision, an Abcam company and was previously called 6527 Protein A/G Magnetic Beads. 6527-1 is the same size as the 1 ml size of ab286842.</p> <p>Protein A/G Magnetic Beads (ab286842) are prepared by covalently coupling Recombinant fusion Protein A/G (contains eight IgG binding domains) to 6% crosslinked magnetically beaded agarose. The coupling technique is optimized to give a high binding capacity for IgG. The capacity of IgG binding is generally greater than 10 mg of human IgG per ml of wet gel.</p> <p>Reusable Protein A/G conjugated magnetic beads with ≥ 10 mg/ml binding capacity for antibody purification from multiple sources.</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.02% Sodium azide Constituent: 99% PBS
Cellular localization	Secreted
Cellular localization	Secreted

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab286842 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IP		Use at an assay dependent concentration.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors