

Product datasheet

Tween 20 (Polyoxyethelenesorbitan Monolaurate) ab128987

1 References

Overview

Product name	Tween 20 (Polyoxyethelenesorbitan Monolaurate)
Tested applications	Suitable for: IHC-P
General notes	<p>Tween 20 detergent (surfactant) ab128987 formulated to be added to buffers and reagents for immunohistochemistry to decrease background staining and enhance reagent spreading during IHC in automated and manual procedures.</p> <p>Procedure:</p> <p>Suggested to use at concentrations between 0.05% and 0.1% in buffers and reagents.</p> <p><u>IHC Protocol Summary</u></p>

Properties

Form	Liquid
Storage instructions	Store at +4°C.
Storage buffer	Constituent: 100% Tween

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab128987 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors