

## Product datasheet

# Donkey Anti-Chicken IgY H&L (FITC) ab63507

★★★★★ [2 Abreviews](#) [15 References](#) [2 Images](#)

### Overview

<b>Product name</b>	Donkey Anti-Chicken IgY H&L (FITC)
<b>Host species</b>	Donkey
<b>Target species</b>	Chicken
<b>Tested applications</b>	<b>Suitable for:</b> IHC-Fr
<b>Immunogen</b>	Purified chicken IgY isolated from normal chicken egg yolks.
<b>Conjugation</b>	FITC. Ex: 493nm, Em: 528nm

### Properties


<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 7.20 Preservative: 0.01% Thimerosal (merthiolate) Constituents: PBS, 50% Glycerol (glycerin, glycerine)
<b>Purity</b>	Protein L purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab63507 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-Fr	★★★★★ (1)	Use at an assay dependent dilution.

### Images




Immunohistochemistry (Frozen sections) - Donkey Anti-Chicken IgY H&L (FITC) (ab63507)

This image is courtesy of an anonymous Abreview.

**ab13970 staining GFP in E10.5 mouse embryo sections by Immunohistochemistry (IHC-Fr - frozen sections).**

Tissue was fixed with paraformaldehyde, permeabilized with PBS/0.5% Triton-X100 and blocked with protein block for 30 minutes at 25°C. Samples were incubated with primary antibody (1/500) for 3 hours at 25°C. FITC-conjugated donkey anti-chicken IgY H&L (ab63507) (1/500) was used as the secondary antibody.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Donkey Anti-Chicken IgY H&L (FITC) (ab63507)

ab63507 was used at dilution 1/100 with the primary antibody **ab13970** in IHC - Free Floating. See the review on **ab13970**.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors