

Product datasheet

Donkey Anti-Goat IgG H&L (HRP polymer) ab214881

★★★★★ [1 Abreviews](#) [1 Image](#)

Overview

Product name	Donkey Anti-Goat IgG H&L (HRP polymer)
Host species	Donkey
Target species	Goat
Specificity	Based on IEP, ab214881 reacts with heavy (gamma) chains on goat IgG, light chains on all goat immunoglobulins. No cross-reactivity is observed to non-immunoglobulin goat serum proteins, IgG from human, mouse, rabbit or rat.
Tested applications	Suitable for: Immunoelectrophoresis, IHC-P
Immunogen	Full length protein corresponding to Goat IgG. (Purified Goat IgG, whole molecule).
Conjugation	HRP polymer

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.1% Proclin 300 Constituents: 1% BSA, 98% PBS
Purity	Affinity purified
Purification notes	ab214881 was affinity purified using solid phase Goat IgG. Purity >95% by SDS-PAGE.
Clonality	Polyclonal
Isotype	IgG


Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab214881 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes

Application	Abreviews	Notes
Immunoelectrophoresis		Use at an assay dependent concentration. This product is provided in a pre-diluted, ready-to-use format. Incubate 150 µl per slide for 1 hour.
IHC-P	★★★★★ (1)	Use at an assay dependent concentration. This product is provided in a pre-diluted, ready-to-use format. Incubate 150 µl per slide for 1 hour.

Images



Human tonsil tissue stained with 1µg/ml goat anti-human IgG (H&L) detected with ab214881 Donkey Anti-Goat H&L (HRP polymer).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Donkey Anti-Goat IgG H&L (HRP polymer) (ab214881)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors