

Product datasheet

Donkey Anti-Mouse IgG H&L (18nm Gold) preadsorbed ab105278

1 References

Overview

Product name	Donkey Anti-Mouse IgG H&L (18nm Gold) preadsorbed
Host species	Donkey
Target species	Mouse
Tested applications	Suitable for: ELISA, WB, Electron Microscopy
Minimal cross-reactivity	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Sheep more details
Immunogen	Full length protein. Mouse IgG whole molecule
Conjugation	Gold 18nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Do Not Freeze.
Storage buffer	pH: 8.50 Preservative: 0.04% Sodium azide Constituents: 1.34% BSA, 0.8% Sodium chloride, 0.25% Sodium borate, 0.09% Dibasic monohydrogen sodium phosphate, 89.08% Water
Purity	Protein G purified
Purification notes	Purified by affinity chromatography.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab105278 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/10 - 1/20.
WB		1/10 - 1/20.
Electron Microscopy		1/10 - 1/20.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors