

Product datasheet

Donkey Anti-Mouse IgG H&L (Alkaline Phosphatase) preadsorbed ab7062

Overview

Product name	Donkey Anti-Mouse IgG H&L (Alkaline Phosphatase) preadsorbed
Host species	Donkey
Target species	Mouse
Tested applications	Suitable for: WB, ICC/IF, Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy
Minimal cross-reactivity	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat, Sheep more details
Immunogen	Full length native mouse IgG (purified).
Conjugation	Alkaline Phosphatase

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.1% Sodium azide Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, 1% BSA
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using mouse IgG coupled to agarose beads.
Conjugation notes	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
Clonality	Polyclonal

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab7062 in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		
ICC/IF		
Dot blot		
ELISA		
IHC-P		
IHC-Fr		
Immunomicroscopy		

Application notes

Dot: Use at an assay dependant dilution.

ELISA: 1/500 - 1/2000.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.

IM: Use at an assay dependant dilution.

WB: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors