

Product datasheet

Donkey Anti-Rabbit IgG H&L (Texas Red ®) ab6800

20 References

Overview

Product name	Donkey Anti-Rabbit IgG H&L (Texas Red ®)
Host species	Donkey
Target species	Rabbit
Tested applications	Suitable for: Immunomicroscopy, Flow Cyt, IHC-P, IHC-Fr, ICC/IF, ELISA
Immunogen	Rabbit IgG, whole molecule
Conjugation	Texas Red ®. Ex: 596nm, Em: 620nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Tripotassium orthophosphate
Purity	Affinity purified
Purification notes	Texas Red Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 1.7 moles Texas Red per mole of Goat IgG
Conjugation notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose
Clonality	Polyclonal

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab6800 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Immunomicroscopy		Use at an assay dependent dilution.
Flow Cyt		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		1/1000 - 1/5000.
ELISA		1/10000 - 1/50000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors