

## Product datasheet

# Goat Anti-Chicken IgA H&L (HRP) ab112817

### 1 References

#### Overview

---

<b>Product name</b>	Goat Anti-Chicken IgA H&L (HRP)
<b>Host species</b>	Goat
<b>Target species</b>	Chicken
<b>Specificity</b>	This antibody may cross react with IgA from other species.
<b>Tested applications</b>	<b>Suitable for:</b> WB, ELISA, IHC-P
<b>Conjugation</b>	HRP

#### Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	Antiserum was solid phase adsorbed to ensure class specificity. ab112817 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to horseradish peroxidase (HRP).
<b>Conjugation notes</b>	Molar enzyme/antibody protein ratio is 4:1.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

#### Applications

---

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab112817 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000 - 1/10000. colorimetric detection; 1/1000 - 1/30000 chemiluminescent detection.

Application	Abreviews	Notes
ELISA		1/10000 - 1/100000. as a primary.
IHC-P		1/200 - 1/500.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors