

Product datasheet

Goat Anti-Chicken IgY H&L ab6872

3 References

Overview

Product name	Goat Anti-Chicken IgY H&L
Host species	Goat
Target species	Chicken
Tested applications	Suitable for: WB, ELISA, ICC/IF, IP, Immunodiffusion, Conjugation, IHC-P, IHC-Fr
Immunogen	Chicken IgY whole molecule

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 0.42% Potassium phosphate
Purity	Affinity purified
Conjugation notes	prepared from monospecific antiserum by immunoaffinity chromatography using Chicken IgG coupled to agarose beads
Clonality	Polyclonal

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab6872 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		
ELISA		
ICC/IF		

Application	Abreviews	Notes
IP		
Immunodiffusion		
Conjugation		
IHC-P		
IHC-Fr		

Application notes suitable for immunoprecipitation, immunodiffusion, conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors