

## Product datasheet

# Goat Anti-Chicken IgY H&L ab97131

## 2 References

### Overview

---

|                            |  |
|----------------------------|--|
| <b>Product name</b>        | Goat Anti-Chicken IgY H&L  |
| <b>Host species</b>        | Goat   |
| <b>Target species</b>      | Chicken  |
| <b>Specificity</b>         | By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgY and with light chains common to other Chicken immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. |
| <b>Tested applications</b> | <b>Suitable for:</b> IHC-P, ELISA, WB, ICC   |

### Properties

---

|                             |  |
|-----------------------------|--|
| <b>Form</b>                 | Liquid   |
| <b>Storage instructions</b> | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| <b>Storage buffer</b>       | pH: 6.8<br>Preservative: 0.1% Sodium azide<br>Constituent: PBS   |
| <b>Purity</b>               | Immunogen affinity purified  |
| <b>Purification notes</b>   | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.                                    |
| <b>Clonality</b>            | Polyclonal   |
| <b>Isotype</b>              | IgG  |

### Applications

---

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab97131 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes   |
|-------------|-----------|---|
| IHC-P       |           | 1/200 - 1/2000.                                     |
| ELISA       |           | 1/1000 - 1/30000. (Primary) 1/100 - 1/500 (Coating) |

| Application | Abreviews | Notes  |
|-------------|-----------|--|
| WB          |           | 1/1000 - 1/10000. (Colorimetric) 1/1000 - 1/30000 (Chemiluminescent) |
| ICC         |           | 1/200 - 1/2000.  |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors