

Product datasheet

Goat Anti-Chicken IgY H&L (FITC) ab6873

★★★★★ [1 Abreviews](#) [4 References](#) [2 Images](#)

Overview

Product name	Goat Anti-Chicken IgY H&L (FITC)
Host species	Goat
Target species	Chicken
Tested applications	Suitable for: Immunomicroscopy, Flow Cyt, ICC/IF, IHC-P, IHC-Fr, ELISA
Immunogen	chicken IgY whole molecule
Conjugation	FITC. Ex: 493nm, Em: 528nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 1% BSA, 0.42% Potassium phosphate
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using chicken IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
Conjugation notes	Diaminotriazinylaminofluorescein (DTAF) (Molecular Weight 530 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.3 moles DTAF per mole of Goat IgG
Clonality	Polyclonal


Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab6873 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Immunomicroscopy		Use at an assay dependent concentration.


Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent dilution.
ICC/IF	★★★★★ (1)	1/1000 - 1/5000.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ELISA		1/10000 - 1/50000.

Images



Immunocytochemistry/ Immunofluorescence - Goat
Anti-Chicken IgY H&L (FITC) (ab6873)

Ab6873 was used at dilution 1/500 with the primary antibody [ab39376](#) in ICC/IF. See the review on [ab39376](#).



Immunohistochemistry (Frozen sections) - Goat
Anti-Chicken IgY H&L (FITC) (ab6873)

Ab6873 was used at dilution 1/500 with the primary antibody [ab13991](#) in IHC-Fr. See the review on [ab13991](#).

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors