

Product datasheet

Goat Anti-Donkey IgG H&L (DyLight® 550) ab98821

1 Image

Overview

Product name	Goat Anti-Donkey IgG H&L (DyLight® 550)
Host species	Goat
Target species	Donkey
Tested applications	<b>Suitable for:</b> Flow Cyt, ICC/IF, IHC-P
Conjugation	DyLight® 550. Ex: 562nm, Em: 576nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
Purity	Immunogen affinity purified
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 550.
Clonality	Polyclonal
Isotype	IgG
General notes	DyLight® 550 replaces DyLight® 549.


Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab98821 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		1/50 - 1/200.
ICC/IF		1/50 - 1/500.

Application	Abreviews	Notes
IHC-P		1/50 - 1/500.

## Images



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors