

## Product datasheet

# Goat F(ab')<sub>2</sub> Anti-Mouse IgG - H&L (AP), pre-adsorbed ab98724

### 1 References

#### Overview

---

<b>Product name</b>	Goat F(ab') <sub>2</sub> Anti-Mouse IgG - H&L (AP), pre-adsorbed
<b>Host species</b>	Goat
<b>Target species</b>	Mouse
<b>Tested applications</b>	<b>Suitable for:</b> ICC, IHC-P, ELISA, WB
<b>Minimal cross-reactivity</b>	Human, Rat <a href="#">more details</a>
<b>Immunogen</b>	Other Immunogen Type corresponding to IgG.
<b>Conjugation</b>	Alkaline Phosphatase

#### Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.58% Sodium chloride, 0.0095% Magnesium chloride, 0.00136% Zinc chloride, 1.19% HEPES, 0.2% BSA
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') <sub>2</sub> fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to Alkaline Phosphatase.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

#### Applications

---

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab98724 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC		Use at an assay dependent dilution.
IHC-P		1/100 - 1/1000.
ELISA		1/1000 - 1/10000. (Primary).
WB		1/1000 - 1/10000. (Chemiluminescent detection).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors