

## Product datasheet

# Goat Anti-Human IgA (SureLight® PE) ab99909

## 1 References

### Overview

---

<b>Product name</b>	Goat Anti-Human IgA (SureLight® PE)
<b>Host species</b>	Goat
<b>Target species</b>	Human
<b>Tested applications</b>	<b>Suitable for:</b> Flow Cyt
<b>Conjugation</b>	SureLight® PE

### Properties

---

<b>Form</b>	Lyophilized: Lyophilized material is stable for up to one year. Reconstitute to 1.0 ml with distilled deionized water vortex gently and allow vial to sit on ice for 20 minutes. After product has been reconstituted, product should be stored at 2-8°C in the dark and be used within 3 months. If further dilution of the conjugate is required, use diluted material within one week.
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 0.58% Sodium chloride, 1.64% Sodium phosphate, 1.71% Sucrose
<b>Purity</b>	Immunogen affinity purified
<b>Conjugation notes</b>	SureLight® R-Phycoerythrin Excitation max: 565-498 nm Emission max.: 578 nm F/P 1.3:1
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

---

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab99909 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		1/100 - 1/1000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors