

## Product datasheet

# Goat Anti-Monkey IgG H&L (FITC) ab112766

## 2 References

### Overview

<b>Product name</b>	Goat Anti-Monkey IgG H&L (FITC)
<b>Host species</b>	Goat
<b>Target species</b>	Monkey
<b>Specificity</b>	By immunoelectrophoresis and ELISA ab112766 reacts specifically with monkey IgG and with light chains common to other monkey immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. ab112766 may cross react with IgG from other species.
<b>Immunogen</b>	Other Immunogen Type corresponding to IgG.
<b>Conjugation</b>	FITC. Ex: 493nm, Em: 528nm

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 99% PBS
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	ab112766 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to fluorescein isothiocyanate (FITC).
<b>Conjugation notes</b>	F/P ratio is 4.2.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors