

## Product datasheet

# Goat Anti-Mouse IgG Fc (DyLight® 488) preadsorbed ab98712

[1 References](#) [1 Image](#)

### Overview

|                                 |   |
|---------------------------------|---|
| <b>Product name</b>             | Goat Anti-Mouse IgG Fc (DyLight® 488) preadsorbed |
| <b>Host species</b>             | Goat  |
| <b>Target species</b>           | Mouse   |
| <b>Tested applications</b>      | <b>Suitable for:</b> Flow Cyt, ICC/IF, IHC-P      |
| <b>Minimal cross-reactivity</b> | Human, Rat <a href="#">more details</a>           |
| <b>Conjugation</b>              | DyLight® 488. Ex: 493nm, Em: 518nm                |

### Properties

|                             |  |
|-----------------------------|--|
| <b>Form</b>                 | Liquid   |
| <b>Storage instructions</b> | Shipped at 4°C. Store at +4°C.   |
| <b>Storage buffer</b>       | pH: 6.8<br>Preservative: 0.09% Sodium azide<br>Constituents: PBS, 0.2% BSA   |
| <b>Purity</b>               | Immunogen affinity purified  |
| <b>Purification notes</b>   | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 488. |
| <b>Clonality</b>            | Polyclonal   |
| <b>Isotype</b>              | IgG  |


### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab98712 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes         |
|-------------|-----------|---------------|
| Flow Cyt    |           | 1/50 - 1/200. |

| Application | Abreviews | Notes         |
|-------------|-----------|---------------|
| ICC/IF      |           | 1/50 - 1/500. |
| IHC-P       |           | 1/50 - 1/500. |

## Images



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors