

Product datasheet

Goat Anti-Mouse IgG H&L (Texas Red ®) ab6787

[69 References](#) [1 Image](#)

Overview

Product name	Goat Anti-Mouse IgG H&L (Texas Red ®)
Host species	Goat
Target species	Mouse
Tested applications	Suitable for: ICC/IF, Flow Cyt, IHC-P, IHC-Fr, ELISA
Immunogen	mouse IgG whole molecule
Conjugation	Texas Red ®. Ex: 596nm, Em: 620nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 6.5 Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Potassium phosphate
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
Conjugation notes	Texas Red™ Sulfonil Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 4.9 moles Texas Redä per mole of Rabbit IgG
Clonality	Polyclonal
Isotype	IgG

Applications


The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab6787 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		1/1000 - 1/5000.
Flow Cyt		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ELISA		1/10000 - 1/50000.

Images



Immunocytochemistry/ Immunofluorescence - Goat
Anti-Mouse IgG H&L (Texas Red ®) (ab6787)

Ab6787 was used at dilution 1/400 with the primary antibody **ab7817** in ICC/IF. See the review on **ab7817**.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors