

Product datasheet

Goat Anti-Mouse IgG1 gamma chain (Cy5 ®) preadsorbed ab136127

4 References

Overview

Product name	Goat Anti-Mouse IgG1 gamma chain (Cy5 ®) preadsorbed
Host species	Goat
Target species	Mouse
Tested applications	Suitable for: ELISA, Flow Cyt, WB, ICC/IF
Minimal cross-reactivity	Human more details
Immunogen	Mouse IgG1 paraproteins
Conjugation	Cy5 ®. Ex: 650nm, Em: 667nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Store In the Dark.
Storage buffer	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: 99% PBS
Purity	Immunogen affinity purified
Purification notes	ab136127 was purified by affinity chromatography on mouse IgG1 covalently linked to agarose and cross-absorbed with Mouse IgM, IgG2a, IgG2b, IgG3 and IgA, pooled Human sera and purified Human paraproteins.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab136127 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.
ICC/IF		Use 0.1µg for 10 ⁶ cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors