

Product datasheet

Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed ab47827

[22 References](#) [1 Image](#)

Overview

Product name	Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed
Host species	Goat
Target species	Mouse
Specificity	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, anti-Peroxidase, Mouse IgG and Mouse IgM.
Tested applications	Suitable for: ICC/IF, ELISA, WB, IHC-P, IHC-Fr, Electron Microscopy
Minimal cross-reactivity	Cow, Horse, Human more details
Immunogen	Full length native proteins (purified) (Mouse)
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA IgG and protease free
Purity	Immunogen affinity purified
Purification notes	Prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption (s) to remove any unwanted reactivities.
Clonality	Polyclonal
Isotype	IgG

Applications


The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab47827 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		1/500 - 1/2500.
ELISA		1/5000 - 1/40000.
WB		1/2000 - 1/10000.
IHC-P		1/500 - 1/2500.
IHC-Fr		1/500 - 1/2500.
Electron Microscopy		Use at an assay dependent dilution.

Images



Ab47827 was used at dilution 1/500 with the primary antibody [ab2869](#) in IHC-P. See the review on [ab2869](#).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed (ab47827)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors