

Product datasheet

Goat Anti-Pig IgG H&L (HRP) ab6915

[3 References](#) [1 Image](#)

Overview

Product name	Goat Anti-Pig IgG H&L (HRP)
Host species	Goat
Target species	Pig
Tested applications	Suitable for: ICC/IF, IHC-P, WB, ELISA, Immunomicroscopy, Dot blot, IHC-Fr
Immunogen	Full length native Swine IgG (purified).
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using pig IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
Conjugation notes	Horseradish Peroxidase (HRP)
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab6915 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		

Application	Abreviews	Notes
IHC-P		
WB		
ELISA		
Immunomicroscopy		
Dot blot		
IHC-Fr		

Application notes

Dot: Use at an assay dependant dilution.

ELISA: 1/5000 - 1/20000.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.


IM: Use at an assay dependant dilution.

WB: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Images



Goat Anti-Pig IgG H&L (HRP) (ab6915) at 1/1000 dilution + Pig IgG, 50 ng

ab6915 incubated for 1 hour at room temperature.

Both the heavy and light chains of pig IgG are detected.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors