

Product datasheet

Goat Anti-Rat IgG H&L (DyLight® 594) preadsorbed
ab98422

1 Image

Overview

Product name	Goat Anti-Rat IgG H&L (DyLight® 594) preadsorbed	
Host species	Goat	
Target species	Rat	
Tested applications	Suitable for: Flow Cyt, ICC/IF, IHC-P	
Minimal cross-reactivity	Chicken, Cow, Goat, Horse, Human, Mouse, Rabbit, Sheep	more details
Conjugation	DyLight® 594. Ex: 593nm, Em: 618nm	

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
Purity	Immunogen affinity purified
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 594.
Clonality	Polyclonal
Isotype	IgG


Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab98422 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		1/50 - 1/200.

Application	Abreviews	Notes
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.

Images



Goat Anti-Rat IgG H&L (DyLight® 594) preadsorbed
(ab98422)

Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors