

Product datasheet

Goat Anti-Syrian hamster IgG H&L (APC) preadsorbed ab130806

2 References

Overview

Product name	Goat Anti-Syrian hamster IgG H&L (APC) preadsorbed
Host species	Goat
Target species	Syrian hamster
Specificity	ab130806 reacts with the heavy and light chains of Syrian Hamster IgG. May also react with the light chains of other Hamster immunoglobulins. Goat Anti-Syrian Hamster IgG H&L (Allophycocyanin) preadsorbed shows minimal cross reactivity with Mouse and Rat immunoglobulins
Tested applications	Suitable for: WB, Dot blot, ICC/IF, Flow Cyt, ELISA
Minimal cross-reactivity	Mouse, Rat more details
Immunogen	Full length protein corresponding to IgG. Golden Syrian Hamster IgG.
Conjugation	APC. Ex: 645nm, Em: 660nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Store In the Dark.
Storage buffer	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: PBS, Sucrose
Purity	Immunogen affinity purified
Purification notes	ab130806 was purified by affinity chromatography on Hamster IgG covalently linked to agarose. ab130806 has been cross absorbed with Mouse and Rat immunoglobulins.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab130806 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
ICC/IF		Use 0.1µg for 10 ⁶ cells.
Flow Cyt		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors