

## Product datasheet

### Anti-Human IgG4 pFc' [HP6023] (Biotin) ab99818

[3 References](#) [1 Image](#)

#### Overview

Product name	Anti-Human IgG4 pFc' [HP6023] (Biotin)
Host species	Mouse
Target species	Human
Specificity	ab99818 reacts with the pFc' portion of the heavy chain of Human IgG4.
Tested applications	<b>Suitable for:</b> ELISA, Dot blot, ICC/IF, IHC-Fr, WB
Immunogen	Full length protein corresponding to Human IgG4. Myeloma protein
Conjugation	Biotin

#### Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Store In the Dark.
Storage buffer	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS
Purity	Protein A purified
Clonality	Monoclonal
Clone number	HP6023
Isotype	IgG3
Light chain type	kappa


#### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab99818 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/1000 - 1/2000.

Application	Abreviews	Notes
<b>Dot blot</b>		Use at an assay dependent dilution. for Dot and slot immunoblotting.
<b>ICC/IF</b>		Use at an assay dependent dilution.
<b>IHC-Fr</b>		Use at an assay dependent dilution.
<b>WB</b>		Use at an assay dependent dilution.

## Images



**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors