

# Mouse Anti-goat IgG for IP (HRP) $\alpha$ b157532

## 7 References

### Overview

Product name	Mouse Anti-goat IgG for IP (HRP)
Host species	Mouse
Target species	Goat
Tested applications	<b>Suitable for:</b> WB
Conjugation	HRP

### Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at -20°C.
Storage buffer	Constituents: 0.16% Sodium phosphate, 0.88% Sodium chloride, 0.1% BSA, 50% Glycerol (glycerin, glycerine)
Clonality	Monoclonal
Clone number	eB270
Isotype	IgG
General notes	Number of blots: At least <b>20 blots</b> (based on a 1 $\mu$ l/ml dilution in 5 ml milk).

#### Important protocol notes:

1. The anti-goat VeriBlot for IP secondary antibody (HRP) detects goat IgG antibodies.
2. The anti-goat VeriBlot for IP secondary antibody (HRP) preferentially detects the non-reduced form of goat IgG over the reduced, SDS-denatured forms.
3. IP sample should be completely reduced/denatured before loaded onto a western blot.
4. Milk should be used as the blocking protein for the immunoblot.

#### Western blot and IP resources:

- a) [Western blot a beginner's guide](#)
- b) [IP protocol](#)

### c) IP troubleshooting tips

## Applications

---

### The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab157532 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration. The dilution will depend on the sensitivity of the HRP substrate. The recommended dilution is 1 µl in 1 ml milk. Based on the recommended dilution in 5 ml milk researchers can perform 20 western blots. This product is recommended for the western blot detection of IP samples.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors