

Product datasheet

Mouse Anti-Human IgG4 pFc' [HP6023] (HRP) ab99817

[13 References](#) [1 Image](#)

Overview

Product name	Mouse Anti-Human IgG4 pFc' [HP6023] (HRP)
Host species	Mouse
Target species	Human
Tested applications	Suitable for: ICC, ELISA, Dot blot, WB, IHC-Fr
Immunogen	Full length protein corresponding to Human IgG4. Human IgG4 myeloma protein
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.40 Constituents: 50% Glycerol, 50% PBS
Purity	Protein A purified
Clonality	Monoclonal
Clone number	HP6023
Isotype	IgG3
Light chain type	kappa


Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab99817 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC		Use at an assay dependent dilution.
ELISA		1/1000 - 1/2000.

Application	Abreviews	Notes
Dot blot		Use at an assay dependent dilution. for Dot and slot immunoblotting.
WB		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.

Images



ELISA plate was coated with purified human IgG1, IgG2, IgG3, and IgG4. Immunoglobulins were detected with serially diluted ab99817.

ELISA - Mouse Anti-Human IgG4 pFc' [HP6023]
(HRP) (ab99817)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors