

## Product datasheet

# Mouse monoclonal [4E3] Anti-Human IgG1 hinge heavy chain (HRP) ab99774

★★★★★ [1 Abreviews](#) [17 References](#)

### Overview

|                     |  |
|---------------------|--|
| Product name        | Mouse monoclonal [4E3] Anti-Human IgG1 hinge heavy chain (HRP) |
| Host species        | Mouse  |
| Target species      | Human  |
| Tested applications | <b>Suitable for:</b> ELISA, ICC                                |
| Immunogen           | Other Immunogen Type corresponding to IgG1.                    |
| Conjugation         | HRP  |

### Properties

|                      |   |
|----------------------|---|
| Form                 | Liquid  |
| Storage instructions | Shipped at 4°C. Store at +4°C.                  |
| Storage buffer       | pH: 7.40<br>Constituents: 50% Glycerol, 50% PBS |
| Purity               | Protein G purified                              |
| Clonality            | Monoclonal                                      |
| Clone number         | 4E3   |
| Isotype              | IgG1  |

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab99774 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes                               |
|-------------|-----------|-------------------------------------|
| ELISA       | ★★★★★ (1) | 1/1000 - 1/2000.                    |
| ICC         |           | Use at an assay dependent dilution. |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors