

Mouse monoclonal [KT98] Anti-Rat IgG2b H&L (HRP)
ab106750

2 References 2 Images

Overview

Product name	Mouse monoclonal [KT98] Anti-Rat IgG2b H&L (HRP)
Host species	Mouse
Target species	Rat
Specificity	ab106750 does not cross-react with Rat IgG1, IgG2a, IgG2c and IgM.
Tested applications	Suitable for: WB, ELISA
Immunogen	Rat IgG2b
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.01% Thimerosal (merthiolate) Constituent: PBS
Purity	Protein A purified
Clonality	Monoclonal
Clone number	KT98
Isotype	IgG1

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab106750 in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000.

Application	Abreviews	Notes
ELISA		1/3000.

Images




All lanes : Mouse monoclonal [KT98] Anti-Rat IgG2b H&L (HRP) (ab106750) at 1/1000 dilution

Lane 1 : Rat IgG1

Lane 2 : Rat IgG2a

Lane 3 : Rat IgG2b

Lane 4 : Rat IgG2c



Rat immunoglobulin isotyping by ELISA using Rat primary antibodies and ab106750

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors