

Product datasheet

Mouse monoclonal [SB62 α] Anti-Rabbit IgG light chain ab99695

Overview

Product name	Mouse monoclonal [SB62 α] Anti-Rabbit IgG light chain
Host species	Mouse
Target species	Rabbit
Tested applications	Suitable for: WB, ELISA
Immunogen	Rabbit Light Chain

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 8.20 Constituent: 100% Borate buffered saline
Purity	Protein A purified
Clonality	Monoclonal
Clone number	SB62a
Isotype	IgG2b
Light chain type	lambda

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab99695 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent dilution.
ELISA		Use a concentration of 1 μ g/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors