

## Product datasheet

# Rabbit F(ab')<sub>2</sub> Anti-Dog IgG H&L (Biotin) ab136769

★★★★★ 1 Abreviews

### Overview

Product name	Rabbit F(ab') <sub>2</sub> Anti-Dog IgG H&L (Biotin)
Host species	Rabbit
Target species	Dog
Tested applications	<b>Suitable for:</b> Flow Cyt, IHC-Fr, IHC-P, ICC, WB, ELISA
Immunogen	Highly purified whole immunoglobulin from Dog.
Conjugation	Biotin

### Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Constituents: 50% PBS, 50% Glycerol
Purity	Immunogen affinity purified
Purification notes	ab136769 is purified by affinity chromatography and is thus particularly specific.
Clonality	Polyclonal
Isotype	IgG

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab136769 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt	★★★★★ (1)	Use at an assay dependent concentration.
IHC-Fr		1/50 - 1/200.
IHC-P		1/50 - 1/200.

Application	Abreviews	Notes
ICC		1/50 - 1/200.
WB		1/200 - 1/2000.
ELISA		1/2000 - 1/20000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors