

# Rabbit Anti-Horse IgG H&L ab6700

## 1 References

### Overview

Product name	Rabbit Anti-Horse IgG H&L
Host species	Rabbit
Target species	Horse
Tested applications	<b>Suitable for:</b> ELISA, WB, IHC-P, IHC-Fr, IP, Immunodiffusion, Conjugation, ICC/IF
Immunogen	Horse IgG whole molecule

### Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.01% Sodium azide Constituents: 0.87% Sodium chloride, 0.164% Sodium phosphate
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Horse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
Clonality	Polyclonal

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab6700 in the following tested applications.  
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		
WB		
IHC-P		

Application	Abreviews	Notes
IHC-Fr		
IP		
Immunodiffusion		
Conjugation		
ICC/IF		

#### Application notes

suitable for immunoprecipitation, immunodiffusion, conjugation and most immunological methods requiring lot-to-lot consistency, high titer and specificity

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors