

Product datasheet

Rabbit Anti-Sheep IgG H&L (Alkaline Phosphatase) ab6748

[3 References](#) [1 Image](#)

Overview

Product name	Rabbit Anti-Sheep IgG H&L (Alkaline Phosphatase)
Host species	Rabbit
Target species	Sheep
Tested applications	Suitable for: IHC-P, WB, ELISA, Immunomicroscopy, Dot blot, ICC/IF, IHC-Fr
Immunogen	Sheep IgG whole molecule
Conjugation	Alkaline Phosphatase

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Do Not Freeze.
Storage buffer	pH: 8 Preservative: 0.1% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.788% Tris HCl, 0.0095% Magnesium chloride, 0.001% Zinc chloride, 50% Glycerol (glycerin, glycerine)
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads
Clonality	Polyclonal
General notes	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons).

Applications


The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab6748 in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		
WB		
ELISA		
Immunomicroscopy		
Dot blot		
ICC/IF		
IHC-Fr		

Application notes

suitable for immunoblotting (western or dot blot), ELISA and immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency

Images



Sheep IgG was detected using ab6748 in Dot Blot analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors