

Rabbit Anti-Sheep IgG H&L preadsorbed ab7106

Overview

Product name	Rabbit Anti-Sheep IgG H&L preadsorbed
Host species	Rabbit
Target species	Sheep
Tested applications	Suitable for: WB, ELISA, ICC/IF, IHC-P, IP, Immunodiffusion, Conjugation, IHC-Fr
Minimal cross-reactivity	Human more details
Immunogen	Full length native Sheep IgG (purified).

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.01% Sodium azide Constituents: 0.164% Sodium phosphate, 1.45% Sodium chloride
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose.
Clonality	Polyclonal

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab7106 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		
ELISA		
ICC/IF		

Application	Abreviews	Notes
IHC-P		
IP		
Immunodiffusion		
Conjugation		
IHC-Fr		

Application notes

Conjugation: Use at an assay dependant concentration.

ID: Use at an assay dependant concentration.

IHC-P: Use at an assay dependant concentration.

IHC-Fr: Use at an assay dependant concentration.

IP: Use at an assay dependant concentration.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors