

Product datasheet

Rabbit Anti-Sheep IgM H&L (Alkaline Phosphatase) ab112761

1 References

Overview

Product name	Rabbit Anti-Sheep IgM H&L (Alkaline Phosphatase)
Host species	Rabbit
Target species	Sheep
Specificity	By immunoelectrophoresis and ELISA ab112761 reacts specifically with sheep IgM. May cross react with IgM from other species.
Tested applications	Suitable for: WB, ELISA, IHC-P, ICC/IF
Immunogen	Other Immunogen Type corresponding to IgM.
Conjugation	Alkaline Phosphatase

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 0.58% Sodium chloride, 0.01% Magnesium chloride, 0.001% Zinc chloride, 1.19% HEPES
Purity	Immunogen affinity purified
Purification notes	Antiserum was solid phase adsorbed to ensure class specificity. ab112761 was isolated by affinity chromatography using antigen coupled to agarose beads.
Conjugation notes	Molar enzyme: antibody protein ratio is 1:1.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab112761 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000 - 1/10000. for chemiluminescent detection
ELISA		1/1000 - 1/10000. as a primary
IHC-P		1/100 - 1/1000.
ICC/IF		1/100 - 1/1000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors